RAW-based Workflow:

Unsorted RAW Images --> FastRawViewer --> Adobe Lr/ACR or other RAW Converter



http://www.fastrawviewer.com website: support@fastrawviewer.com email:

Memory Card

open the folder to be sorted in FastRawViewer

Drag-and-drop the folder on FRV icon or into the FRV window

> (for more options see File menu)

FastRawViewer

Initial Culling [1], [2]

Choose "keepers" and copy them to the designated folder or folders on a hard drive

copy chosen images one-by-one to a designated folder on a hard drive

> **C** — new folder **Shift-C** — same folder

"Keepers" Folder on the Hard Drive

open the folder with the "keepers" in FastRawViewer File --> Recent Folders --> top positions in the drop-down list

FastRawViewer

Initial Adjustments [3]

- 1. Orientation
- 2. Exposure (brightness)
- 3. White Balance
- 4. Contrast (Tone) Curve
- 5. Shadow Boost

Assigning Labels and Ratings

Add Titles and Descriptions

User adjustments are recorded as XMP sidecar files in the current folder with images

run external program directly from **FastRawViewer**

R (Alt-R, Ctrl-R)

if using _Rejected subfolder, make sure to uncheck "Include Subfolders" box in Adobe Lr "From" field

External Program [4] (Adobe Lr/ACR, other RAW converter, RawDigger, etc.)

Catalogue Converting Editing Researching etc.

The Source Folder on a Hard Drive

open the source folder in FastRawViewer

Command-Shift-O (for more options

see File menu)

Initial Culling [1]

Cull rejected images and move them into _Rejected Subfolder of the Source Folder

FastRawViewer

move every rejected image to the _Rejected Subfolder

move current image Alt-Shift-Delete

move previous image Shift-Command-[

The Source Folder on the Hard Drive

Rejected SubFolder of the Source Folder)

in FastRawViewer

files from the _Rejected subfolder won't be opened

continue to work with the Original Folder

restore from Rejected Shift-Command-1

later you can recheck rejected files in order to make sure that you really don't want to keep them

Culling Rejected [1]

Check rejected images and if necessary move them from the _Rejected sublfolder back to the Original Folder or to the _Rejected subfolder of the _Rejected subfolder

FastRawViewer

FRV Technical Tools for Culling based on RAW data:

- 1. RAW Histogram
- 2. OE-, UE indication
- 3. Per-channel View
- 4. Focus Peaking

For a quick start, please look at the "Tips" pop-up at every launch of FRV; for detailed instructions, please read the manual.

[2]

Copy or Move?

Instead of copying "keepers" you can move them to the designated folder:

M - new folder

Shift-M - same folder as previous

Initial Adjustments to an Image **Assigning Ratings and Labels**

For a quick start, please look at the "Tips" pop-up at every launch of FRV; for detailed instructions, please read the manual.

Integration with External Programs

Please look at the "Tips" and read the manual to get detailed instructions how to:

- 1. add an external program to FRV;
- 2. launch an external program;
- 3. adjust settings to save changes in XMP sidecar files;
- 4. match an image display.

